



Product Overview

Scientists have discovered that the earth itself exist of a huge electromagnetic field, but unevenly distributed, only under this huge electromagnetic field, gestated the origin of life. So to speak, no earth's electromagnetic field there would be no any life. As we all know, everyone's body carries definite potential, higher number on children body, approximately 6 volts, the middle-aged people approximately 3-4 volts, the elderly generally about 2-3 volts. With the growth of age or illness, the potential on the human body will gradually reduce.

According to this reality, so there exists with ECG, EEG, electrogastrogram and so forth apparatus applied for detecting all kinds of potential condition for the human. If people get sick, his electrograph will also change correspondingly, usually potential reduced or disorder, if the potential carried by the human is zero, then his life will end likewise. Therefore the human body need to recharge, complementing food is an accustomed charging habits, according to the relationship between electric field and organisms ion, U.S. scientists Franklin made a lot of clinical research and invented a high potential therapy. Use of high-potential therapy for charging on the human body is very worthy of promotion, and is as well an old-line and novel therapy.

High electrical potential regimen apparatus is a new multi-functional machine developed by adopting the latest electronic technology, electromagnetic theory, modern physical therapy. Combining high potential, infrared functions in one, carry through comprehensive regulation and balance for the human cells, tissues, organs systems in all directions from inside to outside, make people eliminate pain and maintain youthful vitality, integrate health care and beauty in one, utilizing high potential therapy apparatus, can achieve both perfect effect of leisure and maintenance, showing the new ideas of high-tech, physical-oriented, green and health.

Product Features

1. adopt micro controller to control, stable treatment performance, low failure rate.
2. high-voltage generating device is manufactured by adopting special techniques, which improve the stability and security of the complete machine as well as features with high-voltage alternating electric field, infrared multiple functions.
3. simple manipulation, intuitive, easy to grasp.
4. features with perfect security protection system on the complete machine, the treatment process is safe and comfortable.
5. rational structure of the complete machine with beautiful sculpt.

Technical parameters

Power: 220 V \pm 10%, 50Hz (110V, subject to the nameplate labeled on the equipment)
Power Consumption: = 100 W
Timing: 30 minutes \pm 3 minutes
Warm Roast temperature: high: 45 degrees, medium: 40 degrees, low: 30 degrees, (\pm 3 degrees), optional
Output voltage: high-grade: 15 KV, medium-grade: 12 KV, low-grade: 9 KV, (\pm 20%)
Automatic adjustment for high-grade, medium-grade and low-grade, circulation frequency 30 seconds.
Fuses: RP1-20, 1A/250V, Φ 5 \times 20
Single Size: (38 \times 35 \times 9) cm (long \times width \times height) cm
Single weight: 5 kg (NW)
Treatment side heavy-load: 120 kg
Instrument Category: electrical appliance safety II, B-type equipment, general equipment
Use environment: humidity: 5 degrees to 40 degrees humidity: \leq 80%,
Atmospheric pressure range: 500 ~ 1060 hPa
Storage environment: temperature: 0 degrees to 40 degrees humidity: \leq 80%,
Atmospheric pressure range: 500 ~ 1060 hPa

The main structure, working principle and effect

Configuration: regimen apparatus is mainly composed by host computer, application segment and so forth

Principle: adopting micro-controller to control the input and output, the data of high voltage part is collected by specialized data acquisition system, the electrical circuit part is used to achieve the function of various potential, infrared, and others required by the treatment. High potential therapy apparatus generates safe high-voltage electrostatic field higher than 10,000 volts, less than the hundreds of micro amp features with constant flow between the positron and negatron, after the body enters into this electrostatic field, all the circulatory system on the human body will be driven to circulate. Chinese medicine "Danxi Xinfu," pointed out: "Blood will run when meet hot" and "condensate when meet cold." Applying warm to dispel cold, promoting blood circulation by removing blood stasis is a gem for China's ancient medicine.

Effect: promote metabolism, regulate autonomic nerve and restore the function of brain cell, improve the blood supply for heart-cerebral blood vessel. Promote blood circulation, regulate blood vessel tension, lower blood viscosity, and improve atherosclerosis and so on. If unstuffed there will be no pain, stuffed will lead to pain. All the diseases and troubles caused by stuffed will mark obvious relieving. In particular the basic therapy of foot induction, because the soles of the feet has corresponding reflective band on each part of the body, known as the second human heart. As we all know, "the therapy of cold from feet bottom", using high-temperature warm roast + high potential therapy can attain very effective treatment effect for the therapy on the sole of the foot.

Attention

1. as for initial users, is supposed to deepfreeze use if encounters such as rapid heartbeat or dizziness case.
2. individual local part have the symptoms such as ants crawling feeling, itching, acid, tingle, mild tingle and so forth are all normal. If the case comparatively serious, you should suspend the use.
3. the people who wear heart pacemaker or other electronic monitoring apparatus, the people who have high fever and various hemorrhage diseases, acute infectious diseases, malignancy, serious heart cardiovascular organic diseases, women during pregnancy period, the people who are serious exhaustion in heart and lung and kidney function, the people who just get through heart surgery and so on, are prohibited from using the equipment.
4. Children, the mentally handicapped people, the elderly and unresponsive people are supposed to use under proper care.
5. In the course of use, user's response should be inquired, if feel too hot, it is necessary to appropriately adjust the level of warm-roast signal. Recommendations: diabetic patients should use the middle or low grade of warm-roast signal.

Maintenance and repairs

1. When clean the apparatus appearance, wipes with wringing-dry towels, and use only after thoroughly dry. Do not use alcohol kinds solvent. Specially notice prevent the equipment from immersing into water. Power supply should be unplugged when clean.
2. In case of thunder lightning, must immediately shut down power, suspend use.
3. Handle carefully, if find cracks on the equipment shell, suspend use.
4. If there are faults and please in time contact with sale agent or contact our factory, unauthorized teardown without permission is prohibited.
5. The equipment can not be covered thoroughly when use it in order to avoid burning caused by overheating.
6. Unplug the power cord if not using in a long time.



Healthy Life »»



Buttocks



waist



Back



Head, cervix



Feet, at every pore

Operation directions

Tips:

1. In order to facilitate clean or enhance soft feeling, you can underlay the towel on the face of the instrument, effect is the same.
2. Signal is only for the warm-roast temperature, no matter choose high, medium or low signal, high-voltage electrostatic field is always the same high.

Start up:

1. Connect to power supply; press the power switch, timing show 30-minute.
2. Press the start button, switching indicator lights up, signal indicator lights up, second-flashing indicator lights up.
High-grade: 15 KV, mid-grade: 12 KV, low-grade: 9 KV, ($\pm 20\%$)
High-grade - mid-grade - low-grade - automatic adjustment, circulation frequency is 30 seconds.
3. Choose warm-roast temperature, press the signal button,
High-grade: 45 degrees, mid-grade: 40 degrees, low-grade: 30 degrees, (± 3 degrees)
4. Treatment part (as left figure)
5. After 30 minutes automatically cut-off, second-flashing indicator light off, at the same time buzzer warn up. After 8 seconds, the machines automatically stop working, and then go into alert status.
6. If you need treatment in succession, just click the start button.

Shutdown:

Press the power supply switch, all the indicator lights off. Disconnect the power plug if not using in a long time.

Note:

Testing ball is only for testing, demonstration other than possessing treatment function.

1. When undertakes tests, testing ball will be placed on the two locations slightly behind the center of treatment surface.
2. When undertakes tests, the iron pillar inside the testing ball will flop circularly with rhythm in accordance with seven second interval.
3. When undertakes treatment, the testing balls do not have to be placed on the treatment cover. Whether place on or not won't affect treatment effect.

Healthy Life »»

Security and storage

The apparatus equips with safety protector

1. when output high voltage and infrared, indicator light shows;
2. current limiting device has restricted the current under the range of absolute security;
3. boast over-current protection device, when the current unusual, it will automatically cut off potential output;
4. when finish the treatment, automatically cut off power output while accompanying with voice alarm clue;
5. At the state of start-up and reset, therapy machine will automatically cut off power output.

Storage climate, environmental conditions suitable temperature range: 20 degrees~+55 degrees;

Relative humidity range: = 93%;

Atmospheric pressure range: 50 kPa ~ 106kPa.

Encasement list

User is supposed to check the product according to encasement list, when there are changes, forgive not further notice.

- 1, host computer one
- 2, power cord one
- 3, testing ball one
- 4, user manual one copy
- 5, certificate of approval of product one copy
- 6, Product warranty cards one copy

Announcement

1. The performance of products continuously improves, subject to the change without notice.
2. All rights reserved, without permission, no one shall be reproduced in any form.
3. If encounter problems when use it, please immediately contact with your local authorized dealer, but also can directly contact with our factory.

Common Troubleshooting

- 1, cannot start-up
 - a) blackout, power supply of network blackout;
 - b) plug is not seated properly or damage;
 - c) power fuse (1 A) damage.
- 2, after start-up, no high-voltage output for the high potential
 - a) not setting strength and weakness for high potential;
 - b) not pressing down startup switch;
 - C) high-voltage fuse (0.5) damage.
- 3, after start-up, no temperature changes for the far-infrared
 - A) not setting high and low for infrared;
 - B) not pressing down startup switch;
 - c) Infrared fuse (2 A) damage.
- 4, the display shows clobber
Network power supply has interference, turn off power and restart at the interval of five minutes.

For other fault, please contact with the maintenance department in our factory or the dealer, must not open the case without authorization, otherwise at your own risk.



Healthy Life »»



Foreword

Dear users:

Welcome to use our company's high electrical potential regimen apparatus! In order to give you a general understanding to the apparatus, we've targeted on this circumstance and configured this brochure for you, please read it prior to operating the health apparatus. As for the presentation of the information, we strive to achieve comprehensiveness and easy to understand, so you can better understand the relevant knowledge of the equipment in order to achieve the purpose of safe, effective and economic use of the regimen apparatus.

Healthy Life »

MIRACLE CO.,LTD

Address: Room B601 Hengsheng Building B , No. 205 Huanshi Zhong Road Guangzhou
China /POST NO.:510010

TEL:86-20-86673107 FAX:86-20-86672750

Email: admin@miracle-china.com

Skype: miracle-china

MSN:china.miracle@hotmail.com QQ:286372662

<http://www.miracle-china.cn>

www.miracle-china.cn