

The User Manual

Ion-Cleanse

Before using the machine, be sure to read through this instruction manual.
This product is designed for only use on healthy person.



The User Manual Ion-Cleanse

Preface

The health and life of the people lie on the volume of electricity in body. If a healthy person, he or she should be taken with 80% negative electricity (usually said negative ions) and 20% positive electricity (usually said positive ions). But the human body is like a giant magnet that attracts and holds onto toxins. Heavy metals, partially oxidized fats, disassimilated proteins, cholesterol deposits, uric acid, plaque, lactic acid and the vast array of chemicals that pervade our water, air and food attach themselves ionized form to joints, organs arteries, nerves and tissues. These toxins disrupt normal body functions creating an environment for disease, allergies, and immune system problems.

Ion generating technology has been around for several decades. It is this technology that is the secret to success in the Ion-Cleanse process by using an alternating positive and negative ion charge during the same treatment. This important discovery now makes it possible to remove unwanted toxic material from everyone regardless of the electrical charges the toxin may hold.

Ions are charged atoms that have gained or lost an electron, which causes them to set up a magnetic field capable of attracting and neutralizing oppositely charged particles and pulling them out of the body through a process called osmosis.

Osmosis is a scientific term used to describe the movement of particles through a cellular membrane from a higher concentration to a lower concentration. The lower concentration in this case refers to the ion field set up by the system in the water. Often times it is possible to reduce pain, swelling, edema and other symptoms caused by a long-term life style and a chronic toxic buildup in the body. The long term effectiveness of the detoxification process depends on the other life-enhancing and life style changes that the patient is willing to make.

What is the Ion-Cleanse?

It is an external detoxification method, which removes toxins through your feet. The body detoxifies while you are comfortably seated with your feet placed in a container of warm water. This is a very efficient way to remove toxic accumulations regardless of where they are located in the body.

The Ion-Cleanse Process generates a stream of positive and negative ions (charged atoms), which attract and attach themselves to oppositely charged toxic particles and draws them out of the body through the skin. Each session

takes about 30 minutes.

The array goes into the water with the hands, feet, or body, and the power supply delivers a small amount of direct current into the array, which causes the metals within it, in combination with water, and salt, to generate positively and negatively charged ions.

The work theory of the Ion-Cleanse

The Ion-Cleanse process is an amazing cleansing method for detoxifying your entire body by bathing your feet in warm water while using positive and negative ions from an ion generator. The process is also known as electrolysis. The array element is the actual ion generator device that sits down into the water and creates the ion charge and flow within that solution. It attaches to the power supply, which controls the array element ion flow and exchange.

Ions are charged atoms that have gained or lost an electron, which causes them to set up a electrostatic field capable of attracting and neutralizing oppositely charged particles and so pulling them out of the body through a process called osmosis. Osmosis is a scientific term used to describe the movement of particles through a cellular membrane from a higher concentration to a lower concentration. The lower concentration in this case refers to the ion field set up by the system in the water. Often at times it is possible to reduce pain and swelling, edema and other symptoms caused by a long-term lifestyle and a chronic toxic build up in the body. The long-term effectiveness of the detoxification process also depends on other life enhancing and life style changes that the patient is willing to make.

Operation

packing list:

Name	Main Machine	Array	Power Cable	Basic Manual	Wrist Band	Metal Box
Quantity	1pcs	1pcs	1pcs	1pcs	1pcs	1pcs

The preparation before use:

- 1) Open the box, take out of the main machine, Array and wrist band, put them in good place for use.
- 2) Connect the power line.
- 3) Connect the array and the wrist band with the main machine.
- 4) Put the array in water.
- 5) It don't need muscle test if not special treatment. Just open the main

machine power.

6) Press the button of "OK" to enter the mode choice; press "+" or "-" to choose the suitable mode. Now there are five modes for choice. Normal people usually choose the third mode.

7) Bath your feet in water after finishing choosing the mode, press "ok" to start setting the time through "+" or "-"; usually set 30 minutes original. If no need to change, press "ok" directly.

8) Re-press "OK" to confirm it and press "START" to start working;

Mode choices:

Though there are usually five modes for choices, because we just use it for normal daily health care, as the common user, suggestion is to choose the third mode with (70+/30- OPTION). It doesn't need the muscle test.

a) All positive electricity: for clinic test usage;

b) All negative electricity: just only for clinic test usage;

Not only all positive but also negative electricity, if use them extremely, it will induce the PH of acid and alkali volume uneven.

c) 70+/30-: It is very popular suitable for most of people;

d) 15+/10-/5+:

e) 10-/15+/5-:

After several thousands muscle test, if the reaction from muscle is positive at first, please choose the mode of "d"; if contrarily, please choose "e".

The operating skills

In spite of any countries, there exists the difference of water quality. Only add a few salt, it can promote the ion-cleanse to conduct electricity. The capability of conducting mainly lies on the water electric ability. There exists the capability of conducting more or less in any kinds of waters except distilled water. The capability of conducting of the water determines the ion effects and impacts the function of detoxification. If the panel shows down 1.8 ampere, it can enhance the capability of conducting if adding a few salts into water. But for maintaining the machine life, please avoid adding salts

The most normal volume display is between 12 and 18 voltage, 1.8 and 2.5 ampere. If down 12 voltage, it easily impacts the measuring volume and leads to be heated.

The basic way of muscle test

Some patients have serious healthy diseases or sensitivity to the magnetic field when seeing the doctor or acanthi treatment, so suggestion is to do the muscle test before use.

1) Bend the arm toward inside or unwind it outside;

2) Will fill a strong force when with pressure, if not, please ask for the special doctor;

3) Put the forefinger on the user's nose and teach him or her to force;

4) Try the arm again, you will find the endurance of user turns more and more weak. If the user also can endure the pressure, it shows no electrode effects with him or her so not to use the machine. Passingly, can ask the user if he or she has medicine termly and if with the mineral and water even;

5) If the user with the reaction of electrode, please directly choose it.

The normal test should be guided by the special the scholar but omitted if family use.

Operating instructions

1. Use ordinary tap water. Do not use distilled water. All water quality is different. Only add a little salt at a time. Adding salt will promote the Ion Detoxification System to conduct electricity. The capability of conducting mainly lies in the waters conductivity quality. The capability of the conducting of the water determines the ion effects and impacts the function of detoxification. If the ampere readout is 1.5 amperes or below, add a little salt to enhance the capability of the conducting. But, for maintaining the machine life, please avoid adding salt

2. Adults need 30 minutes per section

3. 8-15 year old children must use with adult supervision for no more than 15 minutes.

4. Under 50 years, OK to use it every other day for 14 continuous treatments, then stop for three weeks before starting the 14 treatment cycle. Over 50 years, OK to use it every three days for 14 continuous treatments, and then stop for four weeks before starting the 14 treatment cycle. Do not use more than three days in a row. Do not use more than twice per day.

5. The array is needs to be cleaned after every use by soaking in a combination of 30 parts water and 1 part of phosphoric acid or household bleach. Soak at least 5 minutes.

Attention

1. Do not touch power while your hand is wet, please pay attention to your personal safety.

2. Avoid striking the machine.

3. Use clean cloth to polish the dust on the machine's surface, do not use any chemicals.

4. Be gentle when pressing the keys to protect the machine under normal use.

5. You must prop up the machine with the handle while the machine is operating.

Cautions

Please don't use the ion cleanse if the below.

- 1) Who with pacemaker or other magnetic devices.
- 2) Who with no reaction to the muscle test or react weakly to the both poles when foot bath.
- 3) Who have medicine of heart disease termly.
- 4) Who have ever received the surgical of replanting the heart.
- 5) Who have to receive the medicine or mentality treatment.

Remarks

Please note as following:

It should maintain normal blood pressure when any medicine treatments. It must use before the medicine treatment to bath feet by the ion cleanse to keep a natural blood pressure.

If low blood pressure, please have something first. The machine may decrease blood sugar.

The pregnant or lactation women please don't use it, because it maybe dredge out the toxins of the tissue and grease and block producing milk nutrition that is necessary to the infant body. Although not any prove it with bad impaction, but for safety, we suggest no use.

It can be used to the person who is receiving kidney dialysis, hormone insulin medicament and the heart symptom like congestive and half-baked. The machine can help the body to get rid of the toxin that cannot be cleaned by the kidney and heart. It will not impact the medicine treatment and decrease the index of insulin. Someone who accept the replanting surgical of the metal arthrosis, if feel uncomfortable under the magnetic, please stop the foot bath at once.

After use, please make up nutrition absorption such as vitamin, mineral, fatty acid and etc to promote the body metabolism and eject the toxin by itself.

Common troubles

If over 2.8 Amp, it shows the panel wrong with alarm coming. It means high output and then the power off automatically.

If both the current and the voltage are very low, the main machine without good radiation and the panel heavy loaded.

If the power is on, but no information after 30 seconds, for example of the voltage, ampere which will be usually displayed every three seconds. Please refer to the reasons as following:

- 1) The array not in water;
- 2) The electric conduct of the water is not enough, like the distilled water,

suggestion is to add a little salts and mix round to enhance it. When the plate in water, the normal current is between 1.8 and 2.5 ampere, the voltage is between 12 and 18V.

3) The wire of the electrode is not connected with the main machine.

4) There are damages or troubles with the wire of the plate, plug or connecting status.

5) The main machine troubles.

Course suggestion:

1) All people over 3 years old can use the machine if only the adults accompany with them. But the time of foot bath should be shorter a little than the adults. (30 minutes)

2) The people under 50 years old use it every other day. One treatment including the continuous 14 times, after finishing the first treatment, suggestion is continue the treatment after three weeks.

3) The people over 50 years old use it every three days. One treatment including the continuous 14 times, after finishing the first treatment, suggestion is continue the treatment after two weeks.

Maintenance

1) The outer shell of the plate is made up of the acryl that can be turned smooth through bath in water. Adopt 1cm about water and about one spoon of white acid powder, and bath the plate in water a whole night. Suggestion is to clean it after using 5-10 times.

2) The array is easily expendable, along with the different using frequency and habits, the expendable degree also will be different.

3) Please contact with the manufacturer if any other troubles.

Storage

Clean the electrode array before storage. Please note that it cannot be fully immersed into water for washing (keep the cable end dry). You can use a cleaning aid (the distributors can supply) for periodic cleaning of the array. Please keep the array in a dry area when it is not in use.

Do not submerge the bath in water as this may cause the display and electronic components to be damaged.



Color or Particle	Material or Area of the Body Being Detoxified
Black	Detoxifying from liver
Black Flecks	Heavy metals
Brown	Detoxifying from liver, cellular debris, tobacco
Dark Green	Detoxifying from gallbladder
Orange	Detoxifying from joints
Red Flecks	Blood clot material
White Cheese-Like Particles	Most likely yeast
White Foam	Detoxifying from Lymphatic system
Yellow-Green	Detoxifying from kidneys, bladder, urinary tract, female/prostate area

MIRACLE CO.,LTD

Address: Room B601 Hengsheng Building B , No. 205 Huanshi Zhong
Road Guangzhou China /POST NO.:510010

TEL:86-20-86673107 FAX:86-20-86672750

Email: admin@miracle-china.com

Skype: miracle-china

MSN:china.miracle@hotmail.com QQ:286372662

<http://www.miracle-china.cn>

<http://www.miracle-china.cn>