



**Operating instruction
for dual system ion cleanse with mp3**


Operating instruction for dual system ion cleanse with mp3

This advanced hydrotherapy product adds the function of music. You can enjoy music while cleansing your body. Listen to soft beautiful music to relieve tension and put you in a good mood while you detox and cleanse. Say goodbye to the city's uproar, sit down, wear a comfortable earphone and listen to music while enjoying a comfortable ionic hydrotherapy session. Just download your favorite music into the unit! Whatever, you will like to relax yourself! This product has two complete working systems, (you can use two arrays to provide service for two people at the same time). You get double benefit!

Possible ion cleanse Benefits

Assists in restoring body vitality against environmental toxins by eliminating chemicals taken in from drugs, environmental pollution, and foods (particularly processed foods and meats)

- Positive effect on weight loss (Assists in weight loss)
- Headache Relief - Assists in relieving headache problems.
- Slows down aging and improves body flexibility.
- Rejuvenates and energizes your whole body.
- Enhances Nutrient Absorption
- Increase Your Energy
- Improve ability to cleanse body
- Promote metabolism and blood circulation
- Helps Relieves Pain and Tension
- Improved Sleep effectively
- Immune System Booster (Helps Balance Immune System)
- Safe for Children 8 and up



Button function of control panel:

1. Press at all time to increase the volume
Or press to skip to the next song
2. Start switch A (control for user system A)
3. Mode control for user system A (can select mode 1,2, or 3)
4. Show density of ion (amperage) and working time for user system A
5. Show density of ion (amperage) and working time for user system B
6. Mode control for user system B (can select mode 1,2, or 3)
7. Play or stop (music toggle switch)
8. Start switch B (control for user system B)
9. Press at all time to reduce the volume
Or press to skip to the previous song

(The default music volume is very low, you should increase it first.)

Content:

Name	Main Machine	Array	Power Cable	USB cable	Ear phone/ Head phones	Basic Manual	Wrist band
Quantity	1 pcs	2 pcs	1 pcs	1 pcs	2 pcs	1 pcs	2 pcs



Array



Wristband



Earphones



USB cable



Power line

Attention: If missing parts, please contact us.

Now you can give 2 treatments at the same time using just 1 machine! ...

And it also has onboard mp3 music system!

Modern design great looking machine

- Durable hard plastic casing with high quality finish
- Easy to clean
- 220v and 110v (use worldwide)
- Power supply lead (cable & plug for use in your country)
- Treatment time control
- Ion Intensity control
- Water density control
- 2 ion activator probes
- Onboard mp3 music system (plug and play)
- Twin interior fans
- Fuse protected
- Operating Manual with detox information


220v and 110v. Use anywhere in the world!

Portable and easy to operate.

Precautions:

- Persons with a Pacemaker or any other battery-operated, electrical device
- Organ transplant recipients
- Suffer of epileptic seizures
- Pregnant women and breast-feeding mothers
- Do not soak 'open wounds'
- Hypoglycemic person should take food before using
- This device is not recommended for children under 8 years

Operation



1. Make sure the power cable and array are inserted into the host machine socket. Put the array in the washbasin and wear wristband on your hand. Then turn on the power switch at the back of the machine.

2. Press A or B working system switch (3, 4) to select which user's system to set up first.

3. The ion density in the water has differences because of regional water quality; add a small amount of salt to make the ion density rise to 2.0 (The higher the ion density, the more ions are produced; ion density between 1.6 to 2.0 is recommended). The ion density can't exceed 2.4, or else you will hear an alarm, you must add more water.

4.If you would like to listen to music for relaxation or sound therapy, just insert the earphone (the jack lies on the front of machine), press “ play/Stop” . Notice: this is a toggle switch, which has only one button under the word “stop” .

5.Music downloading: connect main machine with computer (USB jack) by USB cable. You will see removable drive device in the “my computer” window of your PC, simply copy your mp3 files to this drive.

Attention:

1. Do not touch power plug if your hands are wet for safety.
2. Avoid striking the machine.
3. Using a clean cloth to clean the machine's surface, do not use any chemical agents..
4. Pressing the keys gently to protect the machine under normal use.

Troubleshooting:

Malfunction	Technician
LCD screen does not illuminate	Check On/Off power switch, power cord connection, blown fuse
The ion probe does not work	Check the connection, Clean thoroughly, Replace the array
Pressing any key does not produce a response.	Turn off power switch for 10 seconds, and then restart.

Warranty Information Card

This warranty covers defects in material and workmanship for one year from the original date of delivery. We will replace any part which in our opinion is defective; if the housing has not been subjected to tampering, alteration, or improper use after delivery and has not been repaired by anyone except us. Our obligation does not include the cost of transportation. We are not responsible for damage in transit, and claims for such damage should be presented to the carrier by the customer. The housing is not warranted against misuse or use in abnormal temperature conditions.

For servicing under this warranty, return any defective part to YOUR DISTRIBUTOR within the period listed below.

Warranty Registration Card

To register your warranty for this product, complete the information below.

Name

Address

City/State/Zip

Telephone

Date Purchased

Colour of Water & Particulate After Cleansing & Toxin Removal

Water Colour	Toxin accumulated area / conditions
Black	Liver, alcohol, asthma, heavy metals
Grey	Heavy metals
Brown	Liver, tobacco, fat waste, cellular debris
Green	Kidneys, bladder, urinary system, female illness, or prostate
Red	Thrombus accumulation, easily bruised
Bottle Green	Immune ability, digestive system
Orange	Arthritis, rheumatism, gout, (joints)
White with Bubbles	Lymph system, insomnia, skin allergy
White with Particulate	Ferment incomplete, flatulence, candida, headache, dream-disturbed sleep

Colour changes in the water are just a partial indicator of the toxins that are leaving your body. The unseen microscopic effects are much more profound.

Easy to clean

Fuse protected

Twin interior fans

Ion Intensity control

Water density control

2 ion activator probes

Treatment time control

220v and 110v (use worldwide)

Onboard mp3 music system (plug and play)

Operating Manual with detox information

Durable hard plastic casing with high quality finish



MIRACLE CO.,LTD

Address: Room B601 Hengsheng Building B , No. 205 Huanshi Zhong Road Guangzhou
China /POST NO.:510010

TEL:86-20-86673107 FAX:86-20-86672750

Email: admin@miracle-china.com

Skype: [miracle-china](#)

MSN:china.miracle@hotmail.com QQ:286372662

<http://www.miracle-china.cn>

<http://www.miracle-china.cn>